relaxmyway

# REFLEX DISC

FOTMASSAGE FOOT MASSAGER

# REFLEX DISC FOTMASSAGE FOOT MASSAGER

## FOTENS ZONTERAPI

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
## THE FOOT REFLEXOLOGY

Reflexology has the advantage over an all-body massage due to the fact that there are approx. 72,000 nerve endings of different endings. Precisely performed massages can influence the body in a positive way and enhance personal well-being.


Careful massages of the foot reflex zones can stimulate metabolism and blood circulation among other things, stimulate digestion and support the self-healing process of the body. The more we care about the vitality of our feet, the better it is for our body: The feet are to be seen as a microsystem in which the whole organism is reflected and can be treated individually.

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- > If the feet are neglected, this can have an effect on the entire body.
- > Muscular tension of the body can arise as a result of poor posture. This in turn can have negative effects on the spinal column and the nerve tracts and result in other physical complaints.
- > Sleeping problems and cold feet are usually a consequence of circulatory disorders.
- > Excessive sitting and physical inactivity can cause blood congestion, swollen legs and muscle cramps.
- > Provide your feet with sufficient energy so that the veins embedded in the muscles can transport the blood throughout the entire body again.

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
Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.

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Massage and wellness experience without a plug socket.

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Magnetic therapy generate magnetic alternating fields in your body. Blood circulation and cell metabolism can be influenced positively and oxygen content in cell tissues increases.

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### REFLEXOLOGY & MAGNETIC THERAPY

- Suitable for people with a lack of exercise
- Promotes removal of toxins
- Soothes strained nerves and muscles
- Stimulates digestive problems
- Strengthens the immune system
- Portable and handy for home, office and travel
- Environmentally friendly electricity free

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## Thank you for purchasing a Casada product.

These instructions for use contain important safety and care instructions and explain how to use the product. Please read them carefully and keep the instructions in a suitable location for future reference.

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### SECURITY INSTRUCTIONS

- Keep the Reflex Disc away from heat and sun.
- Wipe the Reflex Disc regularly using a dust-free, lightly moistened cloth.
- Avoid contact with sharp objects such as rings, sharp fingernails or others.
- An improper use of the band will result in a loss of the liability coverage.

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**WARNING:** A plastic bag was used to pack this product. Keep this plastic bag away from babies and children to avoid suffocation. Plastic bags are no toys.

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### USAGE

Place the Reflex Disc on an evenly and flat surface.

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Please step on the Reflex Disc; you can move your feet on the mat randomly.

**Suggestion:** 10 minutes each time. This action can be repeated many times.

**Benefit:** Improves blood circulation and promotes "Chi". Regular use promotes the flow of energy.

**Suitable for:** Elderly people, people with a lack of exercise, cold feet and muscle discomfort.

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Please step on the Reflex Disc using your heels and move on the mat to the left or to the right. **Suggestion:** 10 minutes each time. This action can be repeated many times.

**Benefits:** Good for digestion system.

**Suitable for:** Middle-aged persons, working adults, people who are "Chi"- less.

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Please step on the Reflex Disc using your toes, you can move on the mat to the left or to the right. **Suggestion:** 5 minutes each time. This action can be repeated many times.

**Benefits:** Relieves shoulder- and headaches, helps to relieve stress.

**Suitable for:** Housewives, persons with sleep disorders, stressed adults.

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Use the Reflex Disc while working, watching TV or reading. Place both or just one foot gently and enjoy the soothing effect.

**Suggestion:** 20 minutes or more. This action can be repeated many times.

**Benefits:** Improves blood circulation and relieves digestive problems.

**Suitable for:** Working adults in air-conditioned environment, elderly people and those with overweight and lack of exercise.

### TEKNISK DATA TECHNICAL DATA

Mått measure:

Ø 33,5 cm

Vikt weight:

0,45 kg

Material material:

EVA-Kunststoff

Certifikat certificates:

  
**CASADA**  
healthcare